

## Digital High-Precision Analytical Balance

### Principle:

Measures mass of substances with **ultra-high precision** using an electromagnetic or strain-gauge system. Suitable for **microgram-level weighing**.

**Capacity:** 220 g

**Readability:** 0.00001 g (10 µg)

**Pan Type:** Glass or stainless steel

**Display:** Digital LCD

### Standard Operating Procedure (SOP)

#### 1. Pre-Operation

- Place balance on a **stable, vibration-free surface**
- Ensure balance is **level** using built-in bubble level
- Switch ON and allow **warm-up for 30 minutes**
- Calibrate using **internal or standard calibration weights**

#### 2. Sample Weighing

- Open draft shield doors
- Tare the balance with container if required
- Add sample carefully using spatula or micro-scoop
- Avoid touching sample with hands
- Close draft shield doors while measuring
- Record weight after reading stabilizes

#### 3. Completion

- Remove sample carefully
- Clean pan and surrounding area
- Switch OFF balance (if used occasionally)

### Applications

- Weighing small quantities of powders, liquids, and herbal extracts
- Analytical and formulation research
- Pharmaceutical quality control
- Microanalysis and assay preparation

### Precautions

- Avoid vibrations and air drafts
- Do not overload pan
- Use gloves or tweezers to handle samples

- Keep balance and surroundings dust-free
- Calibrate regularly