

FLAME PHOTOMETER

Principle:

Based on the emission of characteristic radiation when metal ions are excited in a flame. The intensity of emitted light is proportional to the concentration of the element.

Elements Detected:

Na⁺, K⁺, Ca²⁺, Li⁺

Flame Source: LPG / Natural gas + Air

Detector: Photocell / Photomultiplier

Filters: Element-specific optical filters

Standard Operating Procedure (SOP)

1. Instrument Start-Up

- Switch ON main power supply
- Open air supply first, then gas supply
- Ignite flame carefully
- Allow **10–15 minutes** for stabilization

2. Preparation of Standards & Samples

- Prepare standard solutions of known concentrations
- Prepare blank using distilled water
- Dilute samples if required

3. Instrument Calibration

- Aspirate distilled water → set zero
- Aspirate highest standard → set maximum reading
- Aspirate intermediate standards to check linearity

4. Sample Analysis

- Aspirate sample solution
- Select appropriate filter (Na / K / Ca / Li)
- Record emission intensity
- Repeat for accuracy

5. Data Interpretation

- Compare sample readings with standards
- Calculate concentration using calibration curve

6. Instrument Shutdown

- Aspirate distilled water to clean burner
- Turn OFF gas supply first
- Turn OFF air supply
- Switch OFF instrument

Applications

- Estimation of sodium, potassium, calcium
- Electrolyte analysis
- Pharmaceutical & herbal drug analysis
- Water and soil analysis
- Clinical & research laboratories

Precautions

- Ensure proper ventilation
- Never ignite flame without air supply
- Avoid salt deposits on burner
- Clean nebulizer regularly